

CALIBRATING A WARMING CABINET

The following instructions detail the procedure for calibrating the temperature setting for a Warming Cabinet.

Instructions

1. Remove the Outer and Inner Top Panels (Fig.1)

2. Locate the DIAGNOSTIC BUTTON(s) on the back of the Header (Fig. 1).

Note: The example shown in Figure 1 is for a Dual Chamber Warmer. There will be a Diagnostic Button for each individual chamber.




Fig. 1: Remove Outer and Inner Top Panels and locate Diagnostic Button(s)

- On the front of the Header, hold the SET button, the RED up arrow button, and the DIAGNOSTIC BUTTON simultaneously for 3 seconds.

(**Note:** It is necessary to hold down the Diagnostic Button throughout the entire procedure).




Fig. 2: Controller Display, Header Front

- Press the SET button 12 times until IN-B appears on the display screen. Adjust the temperature using the UP and DOWN arrow buttons -
Example - If the calibrated meter says 108°F, and the Warmer says 106°F, change the number to 2 on the Warmer.
 If the calibrated meter says 104°F and the Warmer says 106°F, use the down arrow and change the IN-B to -2.
- Hold the SET button for 3 seconds to return to the main screen. The desired temperature setting should now be displayed.
- Release the Diagnostic Button.
- Repeat the process for the other chamber, if desired.
- Re-attach the Inner Top Panel and Outer Top Panel.