


X-WIDE ARMBOARD - INSTRUCTIONS FOR USE

The X-Wide armboard is an accessory for a SKYTRON surgical table that is used to support the patient's arm. Its 12" width makes it ideal for larger patients. It easily attaches and detaches to any table using standard U.S. side rails. With multiple locking positions abduction and adduction can be quickly adjusted.



Specifications

The armboard provides 180° range of motion. Contains no Natural Rubber Latex(NRL). The device life is five years under normal use. It must be stored between -4°F to +86°F [-20°C to 30°C].

| | |
|----------------|-----------------|
| Width | 12½" [305mm] |
| Height | 2" [51mm] |
| Length | 21½" [533mm] |
| Product Weight | 5½ lbs [2.3 Kg] |

Install the Armboard

1. Place the release trigger clamp over the side rail at approximately a 45 degree angle. (Figure 1).




Figure 1. Install the Armboard

2. Rotate the armboard downward until the release trigger snaps on to the side rail.
3. Ensure the armboard is properly secured.

Position the Armboard

1. Gently lift the back edge of the armboard to prevent binding, and slide the armboard to the desired location on the side rail (Figure 2).




Figure 2. Position the Armboard

2. Pull the directional trigger and rotate the armboard to the desired position. Release directional trigger to lock.
3. Install the arm pad (Velcro side down).



WARNING

Hazard resulting from incorrect use. The maximum load placed on the armboard may not exceed the appropriate proportion of a patient weighing 500 lbs [227kg].

Remove the Armboard

1. Rotate the armboard 90 degrees to the side rail (Figure 2).
2. Press the release trigger and lift the armboard from the table.

Cleaning

Always follow current AORN/EORNA Journal Guidelines to ensure proper cleaning and disinfection procedure.

Always follow product cleaning and disinfecting utilization instructions and warnings indicated by cleaning product manufacturer.



WARNING

Personal injury to patient or staff may result from a lack of proper maintenance of this equipment.

The following procedure should be followed when cleaning the X-Wide Armboard between cases or operations:

1. Place the armboard and the table top in a level position.
2. Ensure that all power is removed from the table. Take preventive measures to avoid spraying directly into connectors or electrical receptacles.



WARNING

Always follow OSHA/EASHW blood-borne pathogens standards for protective clothing, including gloves, masks, and eye protection when cleaning the surgical table.

3. Remove major contaminants from the armboard and table with disposable materials, following appropriate biohazard waste disposal procedures.
4. Remove the pad and place it on a flat surface for cleaning.

CAUTION

DO NOT use cleaners containing phenolics.

5. Use a ready-to-use detergent or one diluted per the manufacturer's recommendations. Ensure that the active ingredients of the product are compatible with the materials of the SKYTRON X-Wide Armboard and table.
6. Apply cleaning fluid liberally to top and sides of each pad and wipe with a clean lint-free cloth.
7. Using a clean, damp, lint-free cloth, wipe the pads to remove the cleaning fluid.

8. Using a clean, dry, lint-free cloth, wipe the pads to remove all moisture.
9. Repeat Steps 6 through 8 to clean the bottom of the pads. Allow to dry.

CAUTION

When using spray cleaners, DO NOT spray fluids directly into electrical receptacles or components.

Before replacing pads on the table, make sure the pads and all adjoining surfaces are completely dry. Moisture trapped between the pads and adjoining surfaces may cause distortion of table tops.

10. When the cleaning procedure is complete, replace all pads and accessories as applicable.

Disinfection

The following antiseptic solutions are approved for use on the table:

- Sodium Hypochlorite 6% diluted to 0.1% (halogen containing compound)
- Hypo Alcohol (iodine color removing agent)
- Chlorhexidine (chlorhexidine gluconate 0.5%)
- Benzalkonium chloride (invert soap 10%)
- Povidone iodine
- Ethanol 80%
- Oxydol (hydrogen peroxide)
- Isopropyl alcohol (IPA) 99.5%

Perform the following steps when disinfecting the armboard.

1. Remove the pad from the armboard.
2. Apply a proper quantity of disinfectant on a clean and lint-free cloth, then wipe the top, sides, and bottom of the pads with the cloth.

NOTICE

Use enough disinfectant on the cloth to ensure it does not dry when wiping.

3. Disinfect the top of the armboard and the side rails using the same procedure.
4. Wipe all parts with dry, clean, and lint-free cloth within 15 minutes after disinfecting.

Replacement Parts




Figure 3. Armboard Components

| NO. | PART NO. | DESCRIPTION | QTY, |
|-----|--------------|--|------|
| 1 | 2-060-06-1 | ARMBOARD, TOP ONLY | 1 |
| 2 | 2-010-07-3 | ARMBOARD TRIGGER ROTATOR | 1 |
| 3 | 2-010-07-9 | ARMBOARD SLIDE BLOCK # 2 | 1 |
| 4 | 2-010-07-8 | X-WIDE ARMBOARD INDEXING BAR | 1 |
| 5 | 2-010-07-2 | ARMBOARD SLIDE BLOCK # 1 | 1 |
| 6 | 2-010-07-10 | WASHER, S.S., 5/16" | 1 |
| 7 | 2-010-07-7 | ARMBOARD RAIL LOCK | 1 |
| 8 | 2-010-07-5 | ARMBOARD MOUNTNG BLOCK | 1 |
| 9 | 2-010-07-11 | SCREW, S.S., BUTTON HEAD SOCKET, 1/4-20 X 5/8" | 1 |
| N/S | 2-010-07-6 | ARMBOARD BLOCK HUB | 1 |
| N/S | 2-060-06-2WS | 2" WELDED SEAM PAD | 1 |
| N/S | 2-060-06-3P | 3" PAD | 1 |

N/S = Not Shown

Label



Revision History

| Date | Revision | Revision History |
|------------|----------|---|
| 06/05/2017 | 0 | Initial Release |
| 09/05/2017 | 1 | Final Page" Updated Emergo Europe address |

Distributed by:

US - SKYTRON

5085 Corporate Exchange Blvd. S.E.
Grand Rapids, MI 49512 (616) 656-2900
www.skytron.us

Manufactured by:

SchureMed

452 Randolph Street
Abington, MA 02351
US - UNITED STATES



Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands



5085 Corporate Exchange Blvd. S.E.
Grand Rapids, MI 49512 • 616.656.2900