




5085 Corporate Exchange Blvd. S.E.  
Grand Rapids, MI 49512 • 616.656.2900

## UNIVERSAL SIDE RAIL CLAMP

P.N. 6-040-01

### INSTRUCTIONS FOR USE:

The SKYTRON Universal Side Rail Clamp can easily be attached and removed from a SKYTRON table with a U.S. standard side rail. Align the Universal Side Rail Clamp above the side rail and lower into place. Insert the desired accessory and tighten the handle to lock into place. The SKYTRON Universal Side Rail Clamp is designed to fit Round 5/8" (16mm) accessories.



### Replacement Parts:


Handle P.N. 6-040-01-1

### Typical Dimensions:

Width (W) 3 1/4" (82mm)  
Length w/Handle (L) 4 1/2" (113.7mm)

### Weight:

1.95 lbs. (.88kg)



### Authorized Distributors:

#### U.S. - SKYTRON

5085 Corporate Exchange Blvd. S.E.  
Grand Rapids, MI 49512

#### EUROPE - SKYTRON EUROPE BV

FLORESSTRAAT 52  
8022 AD ZWOLLE  
THE NETHERLANDS

EC

REP



### Manufactured by:

Mizuho IKA Co., Ltd.  
Tokyo, Japan



TEC-D-0055 REV0  
3/11