




5085 Corporate Exchange Blvd. S.E.
Grand Rapids, MI 49512 • 616.656.2900


SOLID LEG MODEL PART NUMBERS

X-RAY TOPS




MODEL	COMPLETE SET	A HEAD	B BACK	C SEAT	D FOOT/LEG
300	5-010-04-B	5-010-02-1B	5-010-04-2B	5-010-04-3B	5-010-04-4B
1201K, 6001K with KIDNEY BRIDGE (S.N. 1990-7-1226 & PRIOR - Older models without 90° up back section)	5-010-06-B	5-010-02-1B	5-010-10-2B	5-010-06-3B	5-010-02-4B
6500	5-010-09-B	5-010-02-1B	5-010-09-2B	5-010-09-3B	5-010-09-4B
3500 (S.N. 1996-XX-066 & PRIOR)	5-010-15-B	5-010-02-1B	5-010-09-2B	5-010-15-3B	5-010-09-4B
6600, 6600B	5-010-18-B	5-010-02-1B	5-010-09-2B	5-010-09-3B	5-010-18-4B
1201, 5001, 6001 w/out KIDNEY BRIDGE	5-010-02-B	5-010-02-1B	5-010-02-2B	5-010-02-3B	5-010-02-4B
6001K with KIDNEY BRIDGE (S.N. 1990-7-1227 & LATER - Newer models that have 90° up back section)	5-010-10-B	5-010-02-1B	5-010-10-2B	5-010-10-3B	5-010-10-4B
6500HD	5-010-14-B	5-010-14-1B	5-010-14-2B	5-010-14-3B	5-010-14-4B
3500, 3500B (S.N. 1996-XX-067 & LATER)	5-010-16-B	5-010-02-1B	5-010-09-2B	5-010-16-3B	5-010-09-4B
1700	5-010-19-B	5-010-02-1B	5-010-17-2B	5-010-17-3B	5-010-17-4B
6700 SERIES	5-010-21-B	5-010-21-1B	5-010-21-2B	5-010-21-3B	5-010-21-4B
3600 SERIES	5-010-22-B	5-010-21-1B	5-010-21-2B	5-010-22-3B	5-010-21-4B
3501B, 3501C, 3502, 3503	5-010-23-B	5-010-23-1B	5-010-09-2B	5-010-16-3B	5-010-23-4B
6300	5-010-27-B	5-010-23-1B	5-010-23-2B	5-010-23-3B	5-010-23-4B
1602	5-010-26-B	N/A	N/A	N/A	N/A

DIGITAL CASSETTE X-RAY TOP SET PART NUMBERS




MODEL	COMPLETE SET	A HEAD	B BACK	C SEAT	D FOOT/LEG
3503	5-010-30-D	5-010-30-1D	5-010-30-2D	5-010-30-3D	5-010-30-4D
3603	5-010-29-D	5-010-29-1D	5-010-29-2D	5-010-29-3D	5-010-29-4D

SPLIT LEG MODEL PART NUMBERS





MODEL	COMPLETE SET	A HEAD	B BACK	C SEAT	D FOOT/LEG (LEFT)	E FOOT/LEG (RIGHT)
1200, 5000, 6000 w/out KIDNEY BRIDGE	5-010-03-B	5-010-02-1B	5-010-03-2B	5-010-03-3B	5-010-03-4B	5-010-03-5B
1200K, 6000K with KIDNEY BRIDGE	5-010-05-B	5-010-02-1B	5-010-05-2B	5-010-05-3B	5-010-03-4B	5-010-03-5B
1602SL	5-010-28-B	N/A	N/A	N/A	N/A	N/A

X-RAY TOP REPLACEMENT PARTS




SHORT 5-010-01-7

LONG 5-010-01-7L


LONG 5-010-01-7H
FOR
6500HD, 3600B &
6700 SERIES
(Larger Diameter)STANDARD
X-RAY TOP
STANDOFFDIGITAL X-RAY TOP
CLEAR SIDE STAND OFF FOR
3503 & 3603

MODEL	HEAD	BACK/SEAT	LEG	SCREW
3503	5-010-29-5D	5-010-29-6D	5-010-29-7D	A022065
3603	5-010-30-5D	5-010-30-6D	5-010-30-7D	A022065
Qty	1	1	1	A/R

SAMPLE X-RAY TOP INSTALLATION



IMPORTANT
Apply a drop of Loc-Tite to all screws during assembly,
To prevent loosening during use.



Labels

SKYTRON® GRAND RAPIDS, MI • 616.656.2900

X-Ray Top REF

CAUTION NON-STERILE Consult Instructions for use SN

Mizuho Corporation
3-30-13 Hongo, Bunkyo-ku
Tokyo 113-0033
JAPAN 2017

SKYTRON® GRAND RAPIDS, MI • 616.656.2900

X-Ray Top REF

CAUTION NON-STERILE Consult Instructions for use SN

Komponents Laminated Products
N57 W39741 Industrial Road
PO Box 162
Oconomowoc, WI 53066 2017

Revision History

Date	Revision	Revision History
	0	Initial release.
7/16/2015	1	1602 models added, x-ray top installation added
9/29/2015	2	changes to part numbers under model 1700 and 6300 5-010-17-2B, 5-010-17-3B, 5-010-17-3B changed to 5-010-19-2B, 5-010-19-3B, 5-010-19-3B and 5-10-23-2B, 5-010-23-3B changed to 5-10-27-2B, 5-010-27-3B
10/23/2017	3	Revised document layout to current Skytron standard IFU Revised model & part number information to chart format. Page 2. Revised standoffs to include Digital Cassett X-Ray tops. Changed section title X-Ray Top Standoffs w/Screw to X-Ray Top Repladement Parts. Page 3. Added Digital Cassett X-Ray top for models 3503 & 3603. Added page 4 and labels.

Distributed by:

US - SKYTRON
5085 Corporate Exchange Blvd. S.E.
Grand Rapids, MI 49512 (616) 656-2900
www.skytron.us

Manufactured by:

 Mizuho Corporation
3-30-13 Hongo, Bunkyo-ku
Tokyo 113-0033
JAPAN

 **SKYTRON**®
5085 Corporate Exchange Blvd. S.E.
Grand Rapids, MI 49512 • 616.656.2900